

Anti-OtoCadherin-23 CDH23 Antibody

Catalog Number: A02149-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-OtoCadherin-23 CDH23 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-OtoCadherin-23 CDH23 Antibody catalog # A02149-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H251


Technical Details

Immunogen	Recombinant full length Human CDH23.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200

Anti-OtoCadherin-23 CDH23 Antibody (A02149-1) Images



Western blot analysis of Cadherin 23 expression in HEK293T (A), A549 (B), PC12 (C) whole cell lysates.



Immunohistochemistry (IHC) analysis of CDH23 polyclonal antibody

Anti-OtoCadherin-23 CDH23 Antibody

For Research Use Only. Not for use in diagnostic procedures.