

Anti-Sentrin-specific protease 1 SENP1 Antibody

Catalog Number: A02156-1


Overview

| | |
|----------------------|---|
| Product Name | Anti-Sentrin-specific protease 1 SENP1 Antibody |
| Reactive Species | Human, Mouse |
| Description | Boster Bio Anti-Sentrin-specific protease 1 SENP1 Antibody catalog # A02156-1. Tested in WB applications. This antibody reacts with Human, Mouse. |
| Application | WB |
| Clonality | Monoclonal 5B6 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q9P0U3 |

Technical Details

| | |
|----------------------------|---|
| Immunogen | Synthesized peptide derived from the N-terminal region of human SENP1. at AA range: 1-80 |
| Predicted Reactive Species | Bovine, Canine, Equine, Pig, Rabbit |
| Isotype | IgG |
| Form | Liquid |
| Concentration | This antibody's concentration is >0.5mg/ml. |
| Purification | Immunogen affinity purified |
| Suggested Dilutions | <p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB, 1:500-1:2000</p> |

Anti-Sentrin-specific protease 1 SENP1 Antibody (A02156-1) Images



Western Blot (WB) analysis of specific cells using SENP1 Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Sentrin-specific protease 1 SENP1 Antibody