

## Anti-Phospho-Rad17 (S646) Antibody

Catalog Number: A02159S646

### Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Phospho-Rad17 (S646) Antibody                                                                                                                                                                                                                                                                                                                                                                               |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                |
| Description          | Boster Bio Anti-Phospho-Rad17 (S646) Antibody catalog # A02159S646. Tested in WB, ELISA, IHC applications. This antibody reacts with Human, Mouse, Rat.                                                                                                                                                                                                                                                          |
| Application          | ELISA, IHC, WB                                                                                                                                                                                                                                                                                                                                                                                                   |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                       |
| Formulation          | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                       |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                           |
| Uniprot ID           | O75943                                                                                                                                                                                                                                                                                                                                                                                                           |

### Technical Details

|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Immunogen           | Synthesized phospho-peptide around the phosphorylation site of human Rad17 (phospho Ser646)                           |
| Isotype             | IgG                                                                                                                   |
| Form                | Liquid                                                                                                                |
| Concentration       | 1 mg/ml                                                                                                               |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Suggested Dilutions | WB 1:500-2000<br>IHC 1:50-300<br>ELISA 1:2000-20000                                                                   |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Phospho-Rad17 (S646) Antibody

For Research Use Only. Not for use in diagnostic procedures.