

## Anti-HDAC 9 (P1047) Antibody

Catalog Number: A02177


### Overview

Product Name	Anti-HDAC 9 (P1047) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HDAC 9 (P1047) Antibody catalog # A02177. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UKV0


### Technical Details

Immunogen	Synthetic peptide, corresponding to the C-terminal of Human HDAC 9.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

## Anti-HDAC 9 (P1047) Antibody (A02177) Images



Western blot (WB) analysis of HDAC 9 (P1047) polyclonal antibody at 1:500 dilution  
Lane1:Rat Brain tissue lysate  
Lane2:Mouse Brain tissue lysate  
Lane3:SGC7901 cell lysate  
Lane3:HEK293T cell lysate  
Lane3:MCF-7 cell lysate



Immunohistochemistry (IHC) analyzes of HDAC 9 (P1047) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HDAC 9 (P1047) Antibody

For Research Use Only. Not for use in diagnostic procedures.