

Anti-Uteroglobin/CC10/SCGB1A1 Antibody

Catalog Number: A02201-3

Overview

Product Name	Anti-Uteroglobin/CC10/SCGB1A1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Uteroglobin/CC10/SCGB1A1 Antibody catalog # A02201-3. Tested in IHC, ICC, IF, IF applications. This antibody reacts with Human, Mouse.
Application	IF, IHC, ICC
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P11684

Technical Details

Immunogen	E.coli-derived human SCGB1A1 recombinant protein (Position: M1-C90).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 Immunofluorescence, 1:50-400

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Uteroglobin/CC10/SCGB1A1 Antibody

For Research Use Only. Not for use in diagnostic procedures.