

Anti-IL7R Antibody

Catalog Number: A02222-2


Overview

Product Name	Anti-IL7R Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IL7R Antibody catalog # A02222-2. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P16871

Technical Details

Immunogen	Recombinant fusion protein of human CD127/IL7R(NP_002176.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-IL7R Antibody (A02222-2) Images



Western blot analysis of extracts of various cell lines, using CD127/IL7R antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 180s.

1 Publications Citing This Product

1. PubMed ID: 33344233, He K, Liu S, Xia Y, Xu J, Liu F, Xiao J, Li Y, Ding Q, Lu L, Xiang G, Zhan M. CXCL12 and IL7R as Novel Therapeutic Targets for Liver Hepatocellular Carcinoma Are Correlated With Somatic Mutations and the Tumor Immunological Microenvironment. *Front Oncol.* 2020 Dec 4;10

Visit bosterbio.com/anti-il7r-antibody-a02222-2-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL7R Antibody

For Research Use Only. Not for use in diagnostic procedures.