

Anti-Smad6 Antibody

Catalog Number: A02227


Overview

Product Name	Anti-Smad6 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Smad6 Antibody catalog # A02227. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43541


Technical Details

Immunogen	Recombinant fusion protein of human Smad6(NP_005576.3).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:200

Anti-Smad6 Antibody (A02227) Images



Western blot analysis of extracts of various cell lines, using Smad6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.



Immunofluorescence analysis of NIH-3T3 cells using Smad6 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

1 Publications Citing This Product

1. PubMed ID: 27019660, Effect of Kujie Granule on the Expression of TGF- β /Smads Signaling Pathway in Patients with Ulcerative Colitis

Visit bosterbio.com/anti-smad6-antibody-a02227-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Smad6 Antibody

For Research Use Only. Not for use in diagnostic procedures.