

Anti-Hemopexin HPX Antibody

Catalog Number: A02237


Overview

| | |
|----------------------|--|
| Product Name | Anti-Hemopexin HPX Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Hemopexin HPX Antibody catalog # A02237. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat. |
| Application | IF, IHC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P02790 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Recombinant full length Human HPX. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:2000 IHC: 1:50-1:200 IF: 1:10-1:100 |

Anti-Hemopexin HPX Antibody (A02237) Images



Western blot analysis of HPX polyclonal antibody

1 Publications Citing This Product

1. PubMed ID: -, Wu T,Wang X,Ren K,Huang X,Liu J:The Alterations in Methylene Blue/Light-Treated Frozen Plasma Proteins Revealed by Proteomics.Transfus Med Hemother 2021.doi:10.1159/000515119.

Visit bosterbio.com/anti-hemopexin-hpx-antibody-a02237-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Hemopexin HPX Antibody

For Research Use Only. Not for use in diagnostic procedures.