

Anti-SFPQ Antibody

Catalog Number: A02243-1


Overview

Product Name	Anti-SFPQ Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SFPQ Antibody catalog # A02243-1. Tested in WB,IHC,ICC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P23246


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 ICC: 1:10-1:50

Anti-SFPQ Antibody (A02243-1) Images



Western blot analysis of SFPQ on different lysates using anti-SFPQ antibody at 1/1,000 dilution. Positive control: Lane 1: K562 Lane 2: Daudi Lane 3: Mouse testis Lane 4: Mouse lung



Immunohistochemical analysis of paraffin-embedded rat brain tissue using anti-SFPQ antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SFPQ Antibody

For Research Use Only. Not for use in diagnostic procedures.