

Anti-Adenosine receptor A2a ADORA2A Antibody

Catalog Number: A02268


Overview

Product Name	Anti-Adenosine receptor A2a ADORA2A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Adenosine receptor A2a ADORA2A Antibody catalog # A02268. Tested in WB, ICC applications. This antibody reacts with Human, Mouse, Rat.
Application	ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P29274

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human ADORA2A.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 ICC: 1:50-1:200

Anti-Adenosine receptor A2a ADORA2A Antibody (A02268) Images



Western blot (WB) analysis of ADORA2A polyclonal antibody at 1:500 dilution Lane1:Hela cell lysate Lane2:Raw264.7 cell lysate Lane3:H9C2 cell lysate

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Adenosine receptor A2a ADORA2A Antibody

For Research Use Only. Not for use in diagnostic procedures.