

Anti-COMMD1 Antibody

Catalog Number: A02272


Overview

Product Name	Anti-COMMD1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-COMMD1 Antibody catalog # A02272. Tested in ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, ICC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8N668


Technical Details

Immunogen	Recombinant fusion protein of human COMMD1(NP_689729.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	ICC/IF: 1:50-1:100

Anti-COMMD1 Antibody (A02272) Images



Immunofluorescence analysis of L929 cells using COMMD1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U2OS cells using COMMD1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-COMMD1 Antibody

For Research Use Only. Not for use in diagnostic procedures.