

Anti-COMMD1 Antibody

Catalog Number: A02272-1

About COMMD1

Enables several functions, including phosphatidylinositol-3,4-bisphosphate binding activity; phospholipid binding activity; and protein homodimerization activity. Involved in several processes, including Golgi to plasma membrane transport; negative regulation of protein localization to cell surface; and regulation of cellular protein metabolic process. Acts upstream of or within negative regulation of NF-kappaB transcription factor activity. Located in cytosol; endosome; and nucleoplasm. Part of Cul2-RING ubiquitin ligase complex.

Overview

Product Name	Anti-COMMD1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-COMMD1 Antibody catalog # A02272-1. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8N668

Technical Details

Immunogen	E.coli-derived human COMMD1 recombinant protein (Position: A3-S187).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 Immunoprecipitation, 1:50

ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-COMMD1 Antibody

For Research Use Only. Not for use in diagnostic procedures.