

Anti-FGF9 Antibody Picoband® HRP Conjugated

Catalog Number: A02274-HRP

About FGF9

FGF 9, Fibroblast growth factor 9, is a protein that in humans is encoded by the FGF9 gene. The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. The FGF 9 gene contains 3 exons. By radioactive chromosomal in situ hybridization, the FGF 9 gene is mapped to chromosome 13q11-q12. This protein was isolated as a secreted factor that exhibits a growth-stimulating effect on cultured glial cells. In nervous system, this protein is produced mainly by neurons and may be important for glial cell development. Expression of the mouse homolog of this gene was found to be dependent on Sonic hedgehog (Shh) signaling.

Overview

Product Name	Anti-FGF9 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P31371

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human FGF9, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGF9 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.