

Anti-MKP-1/DUSP1 Antibody Picoband® HRP Conjugated

Catalog Number: A02276-1-HRP

About DUSP1

Dual specificity protein phosphatase 1 is an enzyme that in humans is encoded by the DUSP1 gene. The protein encoded by this gene is a phosphatase with dual specificity for tyrosine and threonine. The encoded protein can dephosphorylate MAP kinase MAPK1/ERK2, which results in its involvement in several cellular processes. This protein appears to play an important role in the human cellular response to environmental stress as well as in the negative regulation of cellular proliferation. Finally, the encoded protein can make some solid tumors resistant to both chemotherapy and radiotherapy, making it a target for cancer therapy.

Overview

Product Name	Anti-MKP-1/DUSP1 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P28562

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human MKP-1/DUSP1, which shares 96.3% amino acid (aa) sequence identity with both mouse and rat MKP-1/DUSP1.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MKP-1/DUSP1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.