

## Anti-Mannose Receptor/MRC1 Picoband® Antibody Biotin Conjugated

Catalog Number: A02285-2-Biotin

### About MRC1

The mannose receptor (Cluster of Differentiation 206, CD206) is a C-type lectin primarily present on the surface of macrophages, immature dendritic cells and liver sinusoidal endothelial cells, but is also expressed on the surface of skin cells such as human dermal fibroblasts and keratinocytes. It is mapped to 10p12.33. The recognition of complex carbohydrate structures on glycoproteins is an important part of several biological processes, including cell-cell recognition, serum glycoprotein turnover, and neutralization of pathogens. The protein encoded by this gene is a type I membrane receptor that mediates the endocytosis of glycoproteins by macrophages. The protein has been shown to bind high-mannose structures on the surface of potentially pathogenic viruses, bacteria, and fungi so that they can be neutralized by phagocytic engulfment.

### Overview

Product Name	Anti-Mannose Receptor/MRC1 Picoband® Antibody Biotin Conjugated
Reactive Species	Human, Monkey, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P22897

### Technical Details

Immunogen	E.coli-derived human Mannose Receptor/MRC1 recombinant protein (Position: D21-A1140).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

## 8 Publications Citing This Product

---

1. PubMed ID: 10.1186/s13287-021-02392-9, IL-1beta pre-stimulation enhances the therapeutic effects of endometrial regenerative cells on experimental colitis

2. PubMed ID: 10.1155/2020/2462039, Protective Effects of Two Safflower Derived Compounds, Kaempferol and Hydroxysafflor Yellow A, on Hyperglycaemic Stress-Induced Podocyte Apoptosis via Modulating of Macrophage M1/M2 Polarization

3. PubMed ID: 10.1155/2020/8263598, Protective Effects of Thalidomide on High-Glucose-Induced Podocyte Injury through In Vitro Modulation of Macrophage M1/M2 Differentiation

Visit [bosterbio.com/anti-mannose-receptor-mrc1-picoband-trade-antibody-a02285-2-boster.html](https://bosterbio.com/anti-mannose-receptor-mrc1-picoband-trade-antibody-a02285-2-boster.html) to see all 8 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mannose Receptor/MRC1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.