

Anti-CD206 MRC1 Antibody

Catalog Number: A02285-3


Overview

Product Name	Anti-CD206 MRC1 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-CD206 MRC1 Antibody catalog # A02285-3. Tested in WB applications. This antibody reacts with Rat,Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P22897

Technical Details

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 610-860 of human Mannose Receptor/CD206 (NP_002429.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:100-1:500

Anti-CD206 MRC1 Antibody (A02285-3) Images



Western blot analysis of HUVEC, using CD206 antibody (A8301) at 1:500 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD206 MRC1 Antibody

For Research Use Only. Not for use in diagnostic procedures.