

## Anti-MSR1 Antibody

Catalog Number: A02349


### Overview

Product Name	Anti-MSR1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MSR1 Antibody catalog # A02349. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P21757

### Technical Details

Immunogen	Recombinant fusion protein of human MSR1(NP_619729.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

## Anti-MSR1 Antibody (A02349) Images



Western blot analysis of extracts of various cell lines, using MSR1 Rabbit pAb at 1:1000 dilution. THP-1 cells were treated by LPS at 37°C for 8 hours and PMA/TPA for 30 minutes after serum-starvation for 16-20 hours. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 1s.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-MSR1 Antibody

For Research Use Only. Not for use in diagnostic procedures.