

Anti-CD204/Msr1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A02349-2-Fluoro550

About Msr1

Macrophage scavenger receptor 1, also known as MSR1, is a protein which in humans is encoded by the MSR1 gene. Enables amyloid-beta binding activity; cargo receptor activity; and low-density lipoprotein particle binding activity. Involved in several processes, including endocytosis; positive regulation of cholesterol storage; and positive regulation of macrophage derived foam cell differentiation. Acts upstream of or within lipoprotein transport. Located in plasma membrane. Is expressed in several structures, including alimentary system; limb; neural tube; reproductive system; and yolk sac. Human ortholog(s) of this gene implicated in Barrett's esophagus; arteriosclerosis; and prostate cancer. Orthologous to human MSR1 (macrophage scavenger receptor 1).

Overview

| | |
|----------------------|--|
| Product Name | Anti-CD204/Msr1 Antibody Picoband® Fluoro550 Conjugated |
| Reactive Species | Mouse |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P30204 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived mouse CD204/Msr1 recombinant protein (Position: M1-D451). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD204/Msr1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.