

Anti-VPS34 (phospho-Ser249) PIK3C3 Antibody

Catalog Number: A02356S249


Overview

Product Name	Anti-VPS34 (phospho-Ser249) PIK3C3 Antibody
Reactive Species	Human
Description	Boster Bio Anti-VPS34 (phospho-Ser249) PIK3C3 Antibody catalog # A02356S249. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NEB9

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 700-800 of Human ULK1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-VPS34 (phospho-Ser249) PIK3C3 Antibody (A02356S249) Images



Western blot (WB) analysis of ULK1 (S757) polyclonal antibody at 1:500 dilution Lane1:L02 whole cell lysate(40ug) Lane2:Panc1 whole cell lysate(40ug) Lane3:C6 whole cell lysate(40ug) Lane4:AML-12 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-VPS34 (phospho-Ser249) PIK3C3 Antibody

For Research Use Only. Not for use in diagnostic procedures.