

## Anti-DNA Ligase I LIG1 Antibody

Catalog Number: A02380-2


### Overview

Product Name	Anti-DNA Ligase I LIG1 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-DNA Ligase I LIG1 Antibody catalog # A02380-2. Tested in WB,IHC applications. This antibody reacts with Human,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P18858


### Technical Details

Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 IHC: 1:50-1:200

## Anti-DNA Ligase I (A02380-2) Images



Western blot analysis of DNA Ligase I on different cell lysates using anti-DNA Ligase I antibody at 1/500 dilution. Positive control: Lane 1: Daudi Lane 2: A431 Lane 3: MCF-7



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-DNA Ligase I antibody. Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DNA Ligase I (A02380-2) Antibody

For Research Use Only. Not for use in diagnostic procedures.