

## Anti-PFKFB3 Antibody

Catalog Number: A02382

### Overview


Product Name	Anti-PFKFB3 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-PFKFB3 Antibody catalog # A02382. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16875

### Technical Details

Immunogen	Recombinant fusion protein of human PFKFB3(NP_004557.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:200

## Anti-PFKFB3 Antibody (A02382) Images

---



Immunofluorescence analysis of U2OS cells using PFKFB3 antibody . Blue: DAPI for nuclear staining.

## 1 Publications Citing This Product

---

1. PubMed ID: 32330381, Zhao Q,Li J,Wu B,Shang Y,Huang X,Dong H,Liu H,Chen W,Gui R,Nie X.Smart Biomimetic Nanocomposites Mediate Mitochondrial Outcome through Aerobic Glycolysis Reprogramming: A Promising Treatment for Lymphoma.ACS Appl Mater Interfaces.2020 May 20;12(20):22687-

Visit [bosterbio.com/anti-pfkfb3-antibody-a02382-boster.html](https://bosterbio.com/anti-pfkfb3-antibody-a02382-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PFKFB3 Antibody

For Research Use Only. Not for use in diagnostic procedures.