

Anti-Ferritin heavy chain FTH1 Antibody

Catalog Number: A02401


Overview

Product Name	Anti-Ferritin heavy chain FTH1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Ferritin heavy chain FTH1 Antibody catalog # A02401. Tested in IHC, ICC/IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P02794


Technical Details

Immunogen	Recombinant fusion protein of human Ferritin Heavy Chain(NP_002023.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200 ICC/IF: 1:50-1:200


Anti-Ferritin heavy chain FTH1 Antibody (A02401) Images




Immunohistochemistry of paraffin-embedded rat stomach using Ferritin Heavy Chain Rabbit pAb at dilution of 1:100 .Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunofluorescence analysis of HeLa cells using Ferritin Heavy Chain Rabbit pAb at dilution of 1:200 . Blue: DAPI for nuclear staining.



Immunohistochemistry of paraffin-embedded human liver cancer using Ferritin Heavy Chain Rabbit pAb at dilution of 1:100 .Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded mouse testis using Ferritin Heavy Chain Rabbit pAb at dilution of 1:100 .Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ferritin heavy chain FTH1 Antibody

For Research Use Only. Not for use in diagnostic procedures.