

Anti-Alpha 1 microglobulin/Ambp Antibody Picoband® PE Conjugated

Catalog Number: A02419-4-PE

About Ambp

Protein AMBP is a protein that in humans is encoded by the AMBP gene. This gene encodes a fusion protein that undergoes proteolytic processing to generate two mature proteins: alpha-1-microglobulin (A1m) is a heme-binding plasma glycoprotein of the lipocalin superfamily of proteins that bind to hydrophobic molecules, whereas bikunin belongs to the superfamily of Kunitz-type protease inhibitors. The transgenic mice specifically lacking bikunin, but not A1m, exhibit female infertility and an increased sensitivity to lung metastasis.

Overview

Product Name	Anti-Alpha 1 microglobulin/Ambp Antibody Picoband® PE Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q07456

Technical Details

Immunogen	E.coli-derived mouse Alpha 1 microglobulin/Ambp recombinant protein (Position: A205-S349).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Alpha 1 microglobulin/Ambp Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.