

## Anti-MSH3 (A80) Antibody

Catalog Number: A02455

### Overview

Product Name	Anti-MSH3 (A80) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MSH3 (A80) Antibody catalog # A02455. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P20585


### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 50-100 of Human MSH3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-MSH3 (A80) Antibody (A02455) Images



Western blot (WB) analysis of MSH3 (A80) pAb at 1:500 dilution Lane1:HEK293T whole cell lysate(20ug) Lane2:MCF-7 whole cell lysate(20ug) Lane3:H9C2 whole cell lysate(20ug) Lane4:CT26 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of MSH3 (A80) pAb in paraffin-embedded liver cancer tissue at 1:100.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MSH3 (A80) Antibody

For Research Use Only. Not for use in diagnostic procedures.