

Anti-PRKCA-binding protein PICK1 Antibody

Catalog Number: A02457-1


Overview

Product Name	Anti-PRKCA-binding protein PICK1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PRKCA-binding protein PICK1 Antibody catalog # A02457-1. Tested in WB,IP,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NRD5


Technical Details

Immunogen	Recombinant protein within N-terminal human PICK1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:2,000-1:5,000 IP: 1:10-1:50 FC: 1:50-1:100

Anti-PRKCA-binding protein PICK1 Antibody (A02457-1) Images



Western blot analysis of PICK1 on SK-Br-3 cell lysate using anti-PICK1 antibody at 1/5,000 dilution.



Flow cytometric analysis of 293T cells with PICK1 antibody at 1/100 dilution (purple) compared with an unlabelled control (cells without incubation with primary antibody; yellow). Alexa Fluor 488-conjugated goat anti-rabbit IgG was used as the secondary antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PRKCA-binding protein PICK1 Antibody

For Research Use Only. Not for use in diagnostic procedures.