

Anti-Syntenin-1 SDCBP Monoclonal Antibody

Catalog Number: A02475-1


Overview

Product Name	Anti-Syntenin-1 SDCBP Monoclonal Antibody
Reactive Species	Dog, Human, Pig, Rat
Description	Boster Bio Anti-Syntenin-1 SDCBP Monoclonal Antibody catalog # A02475-1. Tested in WB applications. This antibody reacts with Human, Rat, Dog, Pig.
Application	WB
Clonality	Monoclonal
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O00560

Technical Details

Immunogen	Purified recombinant human Syntenin-1 protein fragments expressed in E.coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified
Suggested Dilutions	WB 1:1000-1:2000

Anti-Syntenin-1 SDCBP Monoclonal Antibody (A02475-1) Images



Western Blot analysis using Syntenin-1 Monoclonal Antibody against HeLa nuclear extract, HeLa cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Syntenin-1 SDCBP Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.