

Anti-Lysine--tRNA ligase KARS Antibody

Catalog Number: A02479-2


Overview

Product Name	Anti-Lysine--tRNA ligase KARS Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Lysine--tRNA ligase KARS Antibody catalog # A02479-2. Tested in WB,IHC,IP,ICC applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15046


Technical Details

Immunogen	Recombinant protein corresponding to C-terminal human Lysyl tRNA synthetase.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 IP: 1:10-1:50 ICC: 1:50-1:100

Anti-Lysine--tRNA ligase KARS Antibody (A02479-2) Images



Western blot analysis of Lysyl tRNA synthetase on SiHa (1) and MCF-7 (2) cell lysate using anti-Lysyl tRNA synthetase antibody at 1/2,000 dilution.



Immunohistochemical analysis of paraffin-embedded rat testis tissue using anti-Lysyl tRNA synthetase antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Lysine--tRNA ligase KARS Antibody

For Research Use Only. Not for use in diagnostic procedures.