

Anti-CDKAL1 Antibody

Catalog Number: A02489

Overview

Product Name	Anti-CDKAL1 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-CDKAL1 Antibody catalog # A02489. Tested in ELISA, WB applications. This antibody reacts with Human, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q5VV42

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human CDKAL. AA range:10-59
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:40000

Anti-CDKAL1 Antibody (A02489) Images




Figure 1. Western blotting validation for Anti-CDKAL1 Antibody A02489

Western Blot (WB) analysis of COLO205 cells using CDKAL1 polyclonal antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CDKAL1 Antibody