

## Anti-PDE5A Antibody Picoband® PE Conjugated

Catalog Number: A02490-2-PE

### About PDE5A

cGMP-specific phosphodiesterase type 5 is an enzyme from the phosphodiesterase class. It is found in various tissues, most prominently the corpus cavernosum and the retina. It has also been recently discovered to play a vital role in the cardiovascular system. Furthermore, PDE5A also plays a role in signal transduction by regulating the intracellular concentration of cyclic nucleotides. This PDE5A gene is mapped to 4q26.

### Overview

Product Name	Anti-PDE5A Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O76074

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human PDE5A, different from the related mouse sequence by eight amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

---

1. PubMed ID: 10.1016/j.acthis.2018.03.002, An immunohistochemical and ultrastructural analysis of the retina in tadalafil (Cialis) treated rats

Visit [bosterbio.com/anti-pde5a-picoband-trade-antibody-a02490-2-boster.html](https://bosterbio.com/anti-pde5a-picoband-trade-antibody-a02490-2-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PDE5A Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.