

Anti-HSPB8/HSP22 Monoclonal Antibody

Catalog Number: A02492-1


Overview

Product Name	Anti-HSPB8/HSP22 Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HSPB8/HSP22 Monoclonal Antibody catalog # A02492-1. Tested in WB applications. This antibody reacts with Human, Rat, Mouse.
Application	WB
Clonality	Monoclonal 2C3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9UJY1

Technical Details

Immunogen	Recombinant Protein of HSPB8/HSP22
Isotype	IgG1
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified
Suggested Dilutions	WB 1:1000-2000

Anti-HSPB8/HSP22 Monoclonal Antibody (A02492-1) Images



Western blot analysis of 293T with HSPB8/HSP22 Mouse mAb diluted at 1:2, 000.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSPB8/HSP22 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.