

Anti-HSPB8/HSP22 Monoclonal Antibody

Catalog Number: A02492-1


Overview

| | |
|----------------------|---|
| Product Name | Anti-HSPB8/HSP22 Monoclonal Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-HSPB8/HSP22 Monoclonal Antibody catalog # A02492-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat. |
| Application | WB |
| Clonality | Monoclonal 2C3 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Mouse |
| Uniprot ID | Q9UJY1 |

Technical Details

| | |
|----------------------------|--|
| Immunogen | Recombinant Protein |
| Predicted Reactive Species | Bovine, Equine, Guinea Pig |
| Isotype | IgG1 |
| Form | Liquid |
| Concentration | 1 mg/mL. |
| Purification | Protein A purified |
| Suggested Dilutions | <p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB, 1:1000-2000</p> |

Anti-HSPB8/HSP22 Monoclonal Antibody (A02492-1) Images



Western Blot (WB) analysis of 293T with HSPB8/HSP22
Mouse Monoclonal antibody diluted at 1:2000.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSPB8/HSP22 Monoclonal Antibody