

Anti-NTAL/LAT2 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A02496-1-Fluoro550

About LAT2

Linker for activation of T-cells family member 2, also known as WBSCR5 or NTAL is a protein that in humans is encoded by the LAT2 gene. It is mapped to 7q11.23. It has been found that NTAL associated with GRB2, SOS1 and GAB1 after B-cell receptor activation. Phosphorylation of NTAL was strongest in the presence of LCK, ZAP70, SYK, or LYNLAT2. NTAL is a structural and possibly a functional homolog of LAT in non-T cells. It also has been found that BLNK and LAT2 in B cells are a functionally equivalent pair to SLP76 and LAT in T cells. This gene consists of at least 14 exons, and its alternative splicing generates 3 transcript variants, all encoding the same protein.

Overview

Product Name	Anti-NTAL/LAT2 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9GZY6

Technical Details

Immunogen	E.coli-derived human NTAL/LAT2 recombinant protein (Position: R46-S217).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NTAL/LAT2 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.