

Anti-GNAI2 Antibody Picoband® Biotin Conjugated

Catalog Number: A02498-3-Biotin

About GNAI2

Guanine nucleotide-binding protein G(i), alpha-2 subunit is a protein that in humans is encoded by the GNAI2 gene. The protein encoded by this gene is an alpha subunit of guanine nucleotide binding proteins (G proteins). The encoded protein contains the guanine nucleotide binding site and is involved in the hormonal regulation of adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene.

Overview

Product Name	Anti-GNAI2 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P04899

Technical Details

Immunogen	E.coli-derived human GNAI2 recombinant protein (Position: R67-N312). Human GNAI2 shares 97.6% and 98% amino acid (aa) sequence identity with mouse and rat GNAI2, respectively.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GNAI2 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.