

## Anti-GNAI2 Antibody Picoband® PE Conjugated

Catalog Number: A02498-3-PE

### About GNAI2

Guanine nucleotide-binding protein G(i), alpha-2 subunit is a protein that in humans is encoded by the GNAI2 gene. The protein encoded by this gene is an alpha subunit of guanine nucleotide binding proteins (G proteins). The encoded protein contains the guanine nucleotide binding site and is involved in the hormonal regulation of adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene.

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-GNAI2 Antibody Picoband® PE Conjugated   |
| Reactive Species     | Human, Mouse, Rat   |
| Application          | Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality            | Polyclonal  |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .  |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.  |
| Host                 | Rabbit  |
| Uniprot ID           | P04899  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | E.coli-derived human GNAI2 recombinant protein (Position: R67-N312). Human GNAI2 shares 97.6% and 98% amino acid (aa) sequence identity with mouse and rat GNAI2, respectively. |
| Cross Reactivity    | No cross reactivity with other proteins.  |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | PE<br>Excitation Wavelength: 566 nm<br>Emission Wavelength: 574 nm  |
| Suggested Dilutions | Optimal dilutions should be determined by end users.  |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GNAI2 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.