

Anti-S100 S100A1 Antibody

Catalog Number: A02503S100-1


Overview

Product Name	Anti-S100 S100A1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-S100 S100A1 Antibody catalog # A02503S100-1. Tested in WB,IP,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse.
Application	Flow Cytometry, IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P23297


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 FC: 1:50-1:100 IP: 1:10-1:50

Anti-S100 S100A1 Antibody (A02503S100-1) Images



Western blot analysis of S100 on mouse heart cells lysates using anti-S100 antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-S100 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-S100 S100A1 Antibody

For Research Use Only. Not for use in diagnostic procedures.