

Anti-Ox40L TNFSF4 Antibody

Catalog Number: A02554


Overview

Product Name	Anti-Ox40L TNFSF4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Ox40L TNFSF4 Antibody catalog # A02554. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P23510


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human TNFSF4. AA range:31-80
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-1:300 IF 1:100-1:300 ELISA 1:10000

Anti-Ox40L TNFSF4 Antibody (A02554) Images



Western Blot analysis of KB cells using Ox40L Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Western Blot analysis of KB cells using Ox40L Polyclonal Antibody. Secondary antibody was diluted at 1:20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ox40L TNFSF4 Antibody

For Research Use Only. Not for use in diagnostic procedures.