

## Anti-SAA2 Antibody

Catalog Number: A02560-1

### About SAA2

Human SAA proteins are a group of apolipoproteins mainly found in the high density lipoprotein (HDL) portion of plasma. The serum amyloid A (SAA) protein is an acute phase apolipoprotein reactant produced mainly by hepatocytes and under the regulation of inflammatory cytokines (PMID: 12011456). Reactive, secondary amyloidosis is characterized by the extracellular accumulation in various tissues of the SAA2 protein (PMID: 1463770). This antibody can recognize SAA1 and SAA2.

### Overview

Product Name	Anti-SAA2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-SAA2 Antibody catalog # A02560-1. Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	PODJI9

### Technical Details

Immunogen	E.coli-derived human SAA2 recombinant protein (Position: L3-E120).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SAA2 Antibody

For Research Use Only. Not for use in diagnostic procedures.