

Anti-GRP75 (S664) HSPA9 Antibody

Catalog Number: A02561S664


Overview

Product Name	Anti-GRP75 (S664) HSPA9 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GRP75 (S664) HSPA9 Antibody catalog # A02561S664. Tested in WB,IHC,IF,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P38646/Q8N1C8


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids C-terminus of Human GRP75.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC/IF: 1:50-1:200 IP: 1:10-1:100


Anti-GRP75 (S664) HSPA9 Antibody (A02561S664) Images



Western blot (WB) analysis of GRP75 (S664) pAb at 1:5000 dilution Lane1:A549 whole cell lysate(20ug) Lane2:HepG2K562 whole cell lysate(20ug) Lane3:The Brain tissue lysate of Mouse(40ug) Lane4:The Brain tissue lysate of Rat(40ug)



Immunohistochemistry (IHC) analyzes of GRP75 (S664) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.



Immunoprecipitation of A549 whole cell lysate using GRP75 (S664) pAb (Sepharose Bead Conjugate) #BD0047(lane 2 and lane 3) and Nonspecific IgG Control (Sepharose Bead Conjugate)#BD0047 (lane 4 and lane 5) .Lane 1 is 20% input. The western blot was probed using GRP75 (S664) pAb #BS1153. "↑"supernatant; "↓"deposition

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GRP75 (S664) HSPA9 Antibody

For Research Use Only. Not for use in diagnostic procedures.