

## Anti-CD13/ANPEP Picoband® Antibody Fluoro488 Conjugated

Catalog Number: A02591-2-Fluoro488

### About ANPEP

Membrane alanyl aminopeptidase (EC 3.4.11.2) also known as alanyl aminopeptidase (AAP) or aminopeptidase N (AP-N) is an enzyme that in humans is encoded by the ANPEP gene. It is mapped to 15q26.1. Aminopeptidase N is located in the small-intestinal and renal microvillar membrane, and also in other plasma membranes. In the small intestine aminopeptidase N plays a role in the final digestion of peptides generated from hydrolysis of proteins by gastric and pancreatic proteases. Its function in proximal tubular epithelial cells and other cell types is less clear. The large extracellular carboxyterminal domain contains a pentapeptide consensus sequence characteristic of members of the zinc-binding metalloproteinase superfamily. Sequence comparisons with known enzymes of this class showed that CD13 and aminopeptidase N are identical. The latter enzyme was thought to be involved in the metabolism of regulatory peptides by diverse cell types, including small intestinal and renal tubular epithelial cells, macrophages, granulocytes, and synaptic membranes from the CNS. Human aminopeptidase N is a receptor for one strain of human coronavirus that is an important cause of upper respiratory tract infections. Defects in this gene appear to be a cause of various types of leukemia or lymphoma.

### Overview

Product Name	Anti-CD13/ANPEP Picoband® Antibody Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P15144

### Technical Details

Immunogen	E.coli-derived human CD13/ANPEP recombinant protein (Position: D148-S966).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.

Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: -, Yan,Q.,Wang,L.,Song,F.et al.Effects of QDs@Gd3+-NGR on targeted fluorescence-magnetic resonance imaging and inhibition of pancreatic cancer cells. Journal of Materials Research 35, 591–599 (2020). <https://doi.org/10.1557/jmr.2020.53>.

Visit [bosterbio.com/anti-cd13-anpep-picoband-trade-antibody-a02591-2-boster.html](https://bosterbio.com/anti-cd13-anpep-picoband-trade-antibody-a02591-2-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD13/ANPEP Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.