

Anti-JAK3 (I781) Antibody

Catalog Number: A02598-2


Overview

Product Name	Anti-JAK3 (I781) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-JAK3 (I781) Antibody catalog # A02598-2. Tested in WB,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P52333


Technical Details

Immunogen	Synthetic peptide, corresponding to Human JAK3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:2000-1:5000 IF: 1:50-1:200

Anti-JAK3 (I781) Antibody (A02598-2) Images



Western blot (WB) analysis of JAK3 (I781) polyclonal antibody at 1:2000 dilution Lane1:SHSY5Y whole cell lysate(40ug) Lane2:SP2/0 whole cell lysate(40ug) Lane3:C6 whole cell lysate(40ug)



Immunofluorescence analysis of U-87MG cells using JAK3 (I781) antibody at dilution of 1:50.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-JAK3 (I781) Antibody

For Research Use Only. Not for use in diagnostic procedures.