

## Anti-Tyrosine-protein kinase JAK3 JAK3 Antibody

Catalog Number: A02598-4


### Overview

Product Name	Anti-Tyrosine-protein kinase JAK3 JAK3 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Tyrosine-protein kinase JAK3 JAK3 Antibody catalog # A02598-4. Tested in WB applications. This antibody reacts with Human,Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P52333

### Technical Details

Immunogen	Synthetic peptide, corresponding to Human JAK3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000

## Anti-Tyrosine-protein kinase JAK3 JAK3 Antibody (A02598-4) Images



Western blot analysis of JAK3 expression in HeLa (A), Raw264.7 (B) whole cell lysates.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tyrosine-protein kinase JAK3 JAK3 Antibody

For Research Use Only. Not for use in diagnostic procedures.