

## Anti-JNK1/2/3 (T183/Y185) MAPK8 Antibody

Catalog Number: A02608T183Y185


### Overview

Product Name	Anti-JNK1/2/3 (T183/Y185) MAPK8 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-JNK1/2/3 (T183/Y185) MAPK8 Antibody catalog # A02608T183Y185. Tested in WB applications. This antibody reacts with Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P45983/P45984/P53779

### Technical Details

Immunogen	Synthetic peptide, corresponding to Human JNK1/2/3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-JNK1/2/3 (T183/Y185) MAPK8 Antibody (A02608T183Y185) Images



Western blot (WB) analysis of JNK1/2/3 (T183/Y185) polyclonal antibody at 1:500 dilution Lane1:The Brain tissue lysate of Mouse(30ug) Lane2:The Liver tissue lysate of Rat(30ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-JNK1/2/3 (T183/Y185) MAPK8 Antibody

For Research Use Only. Not for use in diagnostic procedures.