

Anti-53BP2/TP53BP2 Antibody Picoband® HRP Conjugated

Catalog Number: A02609-3-HRP

About TP53BP2

tumor protein p53 binding protein 2

Overview

Product Name	Anti-53BP2/TP53BP2 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q13625

Technical Details

Immunogen	E.coli-derived human 53BP2/TP53BP2 recombinant protein (Position: Q97-E1096).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.