

## Anti-53BP2/TP53BP2 Antibody Picoband® PE Conjugated

Catalog Number: A02609-3-PE

### About TP53BP2

tumor protein p53 binding protein 2

### Overview

Product Name	Anti-53BP2/TP53BP2 Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IP, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q13625

### Technical Details

Immunogen	E.coli-derived human 53BP2/TP53BP2 recombinant protein (Position: Q97-E1096).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-53BP2/TP53BP2 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.