

## Anti-IFITM1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A02633-1-Fluoro550

### About IFITM1

Interferon-induced Transmembrane Protein 1 (IFITM1), also called Interferon-induced Protein 17 (IFI17). IFITM1 activity is required for primordial germ cells (PGCs) transit from the mesoderm into the endoderm, and that it appears to act via a repulsive mechanism, such that PGCs avoid Ifitm1-expressing tissues. It is mapped to Chr.11 and belongs to the family of interferon-induced transmembrane proteins (Ifitm/mil/fragilis), which encodes cell surface proteins that may modulate cell adhesion and influence cell differentiation. Interferon-inducible membrane proteins of approximately 17 kDa have been suggested to play a role in the antiproliferative activity of interferons based on their pattern of induction in interferon-sensitive and -resistant cell lines and the ability of a membrane fraction enriched in 17-kDa proteins to inhibit cell growth.

### Overview

Product Name	Anti-IFITM1 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4, 0.02% NaN3.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P13164

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human IFITM1.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 2 Publications Citing This Product

1. PubMed ID: 10.1158/1078-0432.CCR-07-4761, IFN-Induced Transmembrane Protein 1 Promotes Invasion at Early Stage of Head and Neck Cancer Progression

2. PubMed ID: 10.1038/sj.onc.1209807, IFITM1 plays an essential role in the antiproliferative action of interferon- gamma

Visit [bosterbio.com/anti-ifitm1-picoband-trade-antibody-a02633-1-boster.html](http://bosterbio.com/anti-ifitm1-picoband-trade-antibody-a02633-1-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IFITM1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.