

## Anti-CACNA1G Antibody

Catalog Number: A02645


### Overview

Product Name	Anti-CACNA1G Antibody
Reactive Species	Human
Description	Boster Bio Anti-CACNA1G Antibody catalog # A02645. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43497

### Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human CACNA1G.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-CACNA1G Antibody (A02645) Images



Western blot (WB) analysis of CACNA1G pAb at 1:500 dilution Lane1:Myla2059 whole cell lysate(40ug) Lane2:HeLa whole cell lysate(40ug) Lane3:A549 whole cell lysate(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-CACNA1G Antibody

For Research Use Only. Not for use in diagnostic procedures.