

Anti-CXCR7 (A356) ACKR3 Antibody

Catalog Number: A02656


Overview

Product Name	Anti-CXCR7 (A356) ACKR3 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-CXCR7 (A356) ACKR3 Antibody catalog # A02656. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P25106

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 325-375 of Human CXCR7.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-CXCR7 (A356) ACKR3 Antibody (A02656) Images



Western blot (WB) analysis of CXCR7 (A356) pAb at 1:500 dilution Lane1:HEK293T whole cell lysate(40ug) Lane2:Panc1 whole cell lysate(40ug) Lane3:C05-7 whole cell lysate(40ug) Lane4:U-87MG whole cell lysate(40ug) Lane5:K562 whole cell lysate(40ug) Lane6:BV2 whole cell lysate(40ug) Lane7:PC12 whole cell lysate(40ug) Lane8:CT26 whole cell lysate(40ug)

1 Publications Citing This Product

1. PubMed ID: 30636642, Li N,Xu H,Ou Y,Feng Z,Zhang Q,Zhu Q,Cai Z.LPS-induced CXCR7 expression promotes gastric Cancer proliferation and migration via the TLR4/MD-2 pathway.Diagn Pathol.2019 Jan 12;14(1):3.doi:10.1186/s13000-019-0780-x.PMID:30636642;PMCID:PMC6330400.

Visit bosterbio.com/anti-cxcr7-a356-ackr3-antibody-a02656-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CXCR7 (A356) ACKR3 Antibody

For Research Use Only. Not for use in diagnostic procedures.