

Anti-CD3-Epsilon (D71) CD3E Antibody

Catalog Number: A02675-2


Overview

Product Name	Anti-CD3-Epsilon (D71) CD3E Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CD3-Epsilon (D71) CD3E Antibody catalog # A02675-2. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07766

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 31-80 of Human CD3-epsilon.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-CD3-Epsilon (D71) CD3E Antibody (A02675-2) Images



Immunohistochemistry (IHC) analyzes of CD3-epsilon (D71) pAb in paraffin-embedded human brain tissue.

1 Publications Citing This Product

1. PubMed ID: -, Zhao, Q., Shi, M., Yin, C. et al. Dual-Wavelength Photosensitive Nano-in-Micro Scaffold Regulates Innate and Adaptive Immune Responses for Osteogenesis. Nano-Micro Lett. 13, 28 (2021). <https://doi.org/10.1007/s40820-020-00540-z>

Visit bosterbio.com/anti-cd3-epsilon-d71-cd3e-antibody-a02675-2-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD3-Epsilon (D71) CD3E Antibody

For Research Use Only. Not for use in diagnostic procedures.