

## Anti-GluR1 (S857) GRIA1 Antibody

Catalog Number: A02677S857

### Overview


Product Name	Anti-GluR1 (S857) GRIA1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GluR1 (S857) GRIA1 Antibody catalog # A02677S857. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P42261

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 830-880 of Human GluR1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

## Anti-GluR1 (S857) GRIA1 Antibody (A02677S857) Images

---



Immunohistochemistry (IHC) analyzes of GluR1 (S857) pAb in paraffin-embedded human breast carcinoma tissue at 1:50, showing cell membrane staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GluR1 (S857) GRIA1 Antibody

For Research Use Only. Not for use in diagnostic procedures.