

Anti-Cytochrome P450 26B1 CYP26B1 Antibody

Catalog Number: A02691


Overview

Product Name	Anti-Cytochrome P450 26B1 CYP26B1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cytochrome P450 26B1 CYP26B1 Antibody catalog # A02691. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NR63

Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 391-440
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:10000-20000

Anti-Cytochrome P450 26B1 CYP26B1 Antibody (A02691) Images



Western Blot analysis of mouse-brain cells using Antibody diluted at 800. Secondary antibody was diluted at 1:20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytochrome P450 26B1 CYP26B1 Antibody

For Research Use Only. Not for use in diagnostic procedures.